

AirCare Max-2

Portable Compression Therapy

Max-2 portable compression therapy helps increase blood circulation in the lower legs, ankles and feet. It is an effective treatment for pain and swelling caused by conditions associated with poor circulation, such as chronic venous insufficiency and varicose veins.

TGA Approved
Medical Device



- ✓ Improve blood circulation and reduce pain and swelling
- ✓ Heal ulcers and wounds caused by blood pooling
- ✓ Prevent deep vein thrombosis



- ✓ Reduce swelling and inflammation
- ✓ Relieve muscle pain
- ✓ Increase flexibility and range of motion



Portable &
Compact



High Quality &
User-Friendly



MAX-2 Portable Compression Therapy

Max-2 is AirCare high quality portable compression therapy system.

It is an effective treatment for pain and swelling caused by conditions associated with poor circulation, such as chronic venous insufficiency and varicose veins. Max-2 increases blood flow, encouraging your body to deliver more oxygen and nutrients and reduce recovery time, pain and swelling.

Pump Features

- ✓ Ultra-Quite Pump System
- ✓ User-Friendly Smart Operation Panel
- ✓ Six Massage Mode
- ✓ Adjustable pressure comfort level
- ✓ Time indicator
- ✓ Massage position digital display
- ✓ Digital pressure level indicator
- ✓ Quick connect and release plug
- ✓ Pause button
- ✓ Modern Portable Design

Pump Specifications

Voltage:	220-240v / 50Hz
Pressure Range:	20-250mmHg
Pump Weight:	1.5kg
Warranty:	2 Year Limited
SKU:	ACMAX2
Pump Dimensions:	285 x 90 x 90mm



AirCare

Alternating Air Mattresses



PremiumLift

TGA Approved
Medical Device

Address: 1 Force Way Pakenham VIC 3810 Australia

T 1800 771 381

E sales@premiumlift.com.au

www.premiumlift.com.au

All parts and data continually undergo further development and may therefore deviate from the details given.

© **Copyright PremiumLift**

Our Full Terms & Conditions are available by request or can be found on our website. PremiumLift reserves the right to alter or amend this document without notice

Distributor Stamp:

