

# PremiumLift

## Non-Clinical Design

**Promote carer and user wellbeing**



### Unique Hidden Castor Central Locking System Explained

We successfully developed a unique Hidden Castors Central Locking System that revolutionised the home care market and truly enhanced patient's recovery whilst promoting emotional comfort.

Deluxe Floor Bed is fitted with four hidden castors, allowing the user to easily move the bed. The castors can be engaged by pressing and holding the handset Castor Button unit the bed reaches the lowest position.

After the bed is moved or repositioned, press handset Lift Button for two seconds to activate central locking system.



Deluxe Floor Bed can be lowered to just 18cm from the floor level and can elevate to an incredible height of 73cm and meet OH&S height standard requirements.

Visit our website to watch our videos

[www.premiumlift.com.au](http://www.premiumlift.com.au)